

## anti- MRPS25 antibody

### Product Information

Catalog No.:	FNab05356
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	mitochondrial ribosomal protein S25
Synonyms:	RPMS25
Observed MW:	18 kDa
UniprotID :	P82663

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-MRPS25 (IP:FNab05356, 4ug;  
Detection:FNab05356 1:500) with HEK-293 cells  
lysate 1600ug.



HeLa cells were subjected to SDS PAGE followed  
by western blot with FNab05356( MRPS25  
Antibody) at dilution of 1:600