

anti- MS4A6A antibody

Product Information

Catalog No.:	FNab05369
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in signal transduction as a component of a multimeric receptor complex.

Immunogen information

Immunogen:	membrane-spanning 4-domains, subfamily A, member 6A
Synonyms:	4SPAN3, CD20L3, MS4A6
Observed MW:	40-45kd
UniprotID :	Q9H2W1

Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



K-562 cells were subjected to SDS PAGE followed by western blot with FNab05369(MS4A6A antibody) at dilution of 1:400