

anti- MSK1 antibody

Product Information

Catalog No.:	FNab05378
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	ribosomal protein S6 kinase, 90kDa, polypeptide 5
Synonyms:	MSK1, MSPK1, RLPK, RPS6KA5, RSK like protein kinase, RSKL, S6K alpha 5
Observed MW:	80 kDa, 110kDa
UniprotID :	O75582

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, WB
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05378(MSK1 antibody) at dilution of 1:1000