

anti- MUC16,CA125 antibody

Product Information

Catalog No.:	FNab05430
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	0A11
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	mucin 16, cell surface associated
Synonyms:	None
Observed MW:	90-100 kDa
Uniprot ID :	Q8WXI7

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution:IHC:	1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab05430(CA125,MUC16 antibody) at dilution of 1:50