

anti- MX1 antibody

Product Information

Catalog No.:	FNab05450
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Interferon-induced dynamin-like GTPase with antiviral activity against a wide range of RNA viruses and some DNA viruses. Its target viruses include negative-stranded RNA viruses and HBV through binding and inactivation of their ribonucleocapsid. May also antagonize reoviridae and asfarviridae replication. Inhibits thogotovirus(THOV) replication by preventing the nuclear import of viral nucleocapsids. Inhibits La Crosse virus(LACV) replication by sequestering viral nucleoprotein in perinuclear complexes, preventing genome amplification, budding, and egress. Inhibits influenza A virus(IAV) replication by decreasing or delaying NP synthesis and by blocking endocytic traffic of incoming virus particles. Enhances ER stress-mediated cell death after influenza virus infection. May regulate the calcium channel activity of TRPCs.

Immunogen information

Immunogen:	Myxovirus resistance protein 1
Synonyms:	IFI 78K, IFI78, Interferon induced protein p78, MX, MX1, MxA, Myxovirus resistance protein 1
Observed MW:	76 kDa
UniprotID :	P20591

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com


Application

Reactivity: Human, Mouse


Tested Application: ELISA, WB, IHC, IP, IF

Recommended dilution: WB: 1:500-1:1000; IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200

Image:





Immunohistochemistry of paraffin-embedded human lung cancer slide using FNab05450(MX1 Antibody) at dilution of 1:50



Immunofluorescence analysis of A549 cells using FNab05450(MX1 Antibody) at dilution of 1:50

IP result of anti-MX1 (FNab05450 for IP and Detection) with A549 cell.



mouse spleen tissue were subjected to SDS PAGE followed by western blot with FNab05450(MX1 antibody) at dilution of 1:300