

anti- MYH6 antibody

Product Information

Catalog No.:	FNab05477
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Muscle contraction.

Immunogen information

Immunogen:	myosin, heavy chain 6, cardiac muscle, alpha
Synonyms:	MYHCA
Observed MW:	220 kDa
UniprotID :	P13533

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:2000-1:20000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab05477(MYH6 Antibody) at dilution of 1:50



rat heart tissue were subjected to SDS PAGE followed by western blot with FNab05477(MYH6 Antibody) at dilution of 1:5000