

anti- NEDD8 antibody

Product Information

Catalog No.:	FNab05647
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

NEDD8 (neural precursor cell expressed, developmentally down-regulated 8) is a protein-coding gene. Diseases associated with NEDD8 include glomuvenous malformation. GO annotations related to this gene include ubiquitin protein ligase binding.

Immunogen information

Immunogen:	neural precursor cell expressed, developmentally down-regulated 8
Synonyms:	NEDD-8
Observed MW:	12 kDa
UniprotID :	Q15843

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human esophagus using FNab05647(NEDD8 Antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab05647(NEDD8 antibody) at dilution of 1:1000