

anti- NELFB antibody

Product Information

Catalog No.:	FNab05657
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Essential component of the NELF complex, a complex that negatively regulates the elongation of transcription by RNA polymerase II. The NELF complex, which acts via an association with the DSIF complex and causes transcriptional pausing, is counteracted by the P-TEFb kinase complex. May be able to induce chromatin unfolding.

Immunogen information

Immunogen:	cofactor of BRCA1
Synonyms:	COBRA1, KIAA1182
Observed MW:	66 kDa
UniprotID :	Q8WX92

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP, IF
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000; IF: 1:10-1:100


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



IP Result of anti-NELFB (IP:FNab05657, 3ug;
Detection:FNab05657 1:500) with mouse brain
tissue lysate 6000ug.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab05657(NELFB antibody) at dilution of 1:1000

