

NFH antibody

Product Information


Catalog No.:	FNab05686
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1C3
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	neurofilament, heavy polypeptide
Synonyms:	200 kDa neurofilament protein, NEFH, NEH, neurofilament, NF, NF H, NF200, NF-H
Observed MW:	200 kDa
Uniprot ID :	P12036

Application

Reactivity:	Human, Pig, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:1000-1:8000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab05686(NF200antibody) at dilution of 1:50



pig brain tissue were subjected to SDS PAGE followed by western blot with FNab05686(NF200 Antibody) at dilution of 1:4000