

## anti- NFKB2,p52,p100-Specific antibody

### Product Information

Catalog No.:	FNab05707
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	nuclear factor of kappa light polypeptide gene enhancer in B-cells 2(p49/p100)
Synonyms:	Nuclear factor NF-kappa-B p100 subunit,DNA-binding factor KBF2,H2TF1,Lymphocyte translocation chromosome 10 protein,Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2,Oncogene Lym-10,NFKB2
Observed MW:	52/110kd
Uniprot ID :	Q00653

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human testis using FNab05707(NFKB2,p52,p100-Specific antibody) at dilution of 1:50



Hela cells were subjected to SDS PAGE followed by western blot with FNab05707(NFKB2,p52,p100-Specific Antibody) at dilution of 1:1000