

anti- NOP2 antibody

Product Information

Catalog No.:	FNab05794
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

S-adenosyl-L-methionine-dependent methyltransferase that specifically methylates the C(5) position of cytosine 4447 in 28S rRNA(Probable). May play a role in the regulation of the cell cycle and the increased nucleolar activity that is associated with the cell proliferation(Probable).

Immunogen information

Immunogen:	NOP2 nucleolar protein homolog(yeast)
Synonyms:	NOL1, NSUN1
Observed MW:	100-120 kDa
UniprotID :	P46087

Application


Reactivity:	Human, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human lymphoma using FNab05794(NOP2 antibody) at dilution of 1:50



A2780 cells were subjected to SDS PAGE followed by western blot with FNab05794(NOP2 Antibody) at dilution of 1:2000