

anti- NOX2 antibody

Product Information

Catalog No.:	FNab05804
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	cytochrome b-245, beta polypeptide
Synonyms:	CYBB, Cytochrome b 245 heavy chain, NOX2, p91 PHOX
Observed MW:	65 kDa
Uniprot ID :	P04839

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab05804(NOX2 Antibody) at dilution of 1:100



mouse spleen tissue were subjected to SDS PAGE followed by western blot with FNab05804(NOX2 antibody) at dilution of 1:1000