

anti- NPR2 antibody

Product Information

Catalog No.:	FNab05823
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Receptor for the C-type natriuretic peptide NPPC/CNP hormone. Has guanylatecyclase activity upon binding of its ligand. May play a role in the regulation of skeletal growth.

Immunogen information

Immunogen:	natriuretic peptide receptor B/guanylatecyclaseB(atrionatriuretic peptide receptor B)
Synonyms:	ANPRB
Observed MW:	111 kDa, 117 kDa
UniprotID :	P20594

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



HEK293 cells were subjected to SDS PAGE followed by western blot with FNab05823(NPR2 antibody) at dilution of 1:1000