

anti- NSE antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab05867 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Protein A+G purification |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Mouse |
| Clonality: | monoclonal |
| Clone ID: | 0F8 |
| IsoType: | IgG2a |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

None

Immunogen information

| | |
|--------------|--|
| Immunogen: | enolase 2(gamma, neuronal) |
| Synonyms: | enolase 2 (gamma, neuronal), Gamma enolase, Neural enolase, Neuron specific enolase, NSE |
| Observed MW: | 47 kDa |
| UniprotID : | P09104 |

Application


| | |
|-----------------------|-------------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IHC |
| Recommended dilution: | WB: 1:2000-1:20000; IHC: 1:20-1:200 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab05867(ENO2 Antibody) at dilution of 1:200



U-251 cells were subjected to SDS PAGE followed by western blot with FNab05867(ENO2 Antibody) at dilution of 1:4000