

anti- NUB1 antibody

Product Information

Catalog No.:	FNab05887
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Specific down-regulator of the NEDD8 conjugation system. Recruits NEDD8, UBD, and their conjugates to the proteasome for degradation. Isoform 1 promotes the degradation of NEDD8 more efficiently than isoform 2.

Immunogen information

Immunogen:	negative regulator of ubiquitin-like proteins 1
Synonyms:	BS4, NEDD8 ultimate buster 1, NUB1, NUB1L, NYREN18
Observed MW:	69 kDa
Uniprot ID :	Q9Y5A7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



IP Result of anti-NUB1 (IP: FNab05887, 3ug;
Detection: FNab05887 1:600) with HEK-293 cells
lysate 2500ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab05887(NUB1 antibody)
at dilution of 1:600