

anti- NUDT6 antibody

Product Information

Catalog No.:	FNab05906
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May contribute to the regulation of cell proliferation.

Immunogen information

Immunogen:	nudix(nucleoside diphosphate linked moiety X)-type motif 6
Synonyms:	FGF2AS
Observed MW:	72 kDa
UniprotID :	P53370

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275


Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human gliomas using FNab05906(NUDT6 antibody) at dilution of 1:50



IP Result of anti-NUDT6 (IP:FNab05906, 4ug; Detection:FNab05906 1:500) with HepG2 cells lysate 4000ug.



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05906(NUDT6 antibody) at dilution of 1:500