

## anti- Occludin antibody

### Product Information

Catalog No.:	FNab05957
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

This gene encodes an integral membrane protein that is required for cytokine-induced regulation of the tight junction paracellular permeability barrier. Mutations in this gene are thought to be a cause of band-like calcification with simplified gyration and polymicrogyria (BLC-PMG), an autosomal recessive neurologic disorder that is also known as pseudo-TORCH syndrome. Alternative splicing results in multiple transcript variants. A related pseudogene is present 1.5 Mb downstream on the q arm of chromosome 5.

### Immunogen information

Immunogen:	occludin
Synonyms:	occludin, OCLN
Observed MW:	59 kDa
Uniprot ID :	Q16625

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:20-1:200

1

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)


Image:



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab05957( Occludin Antibody) at dilution of 1:500



Immunofluorescent analysis of MCF7 cells, using antibody FNab05957 at 1:50 dilution.



HUVEC cells were subjected to SDS PAGE followed by western blot with FNab05957 (Occludin antibody) at dilution of 1:1500