

anti- ODF3 antibody

Product Information

Catalog No.:	FNab05973
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	outer dense fiber of sperm tails 3
Synonyms:	hSHIPPO, ODF3, Outer dense fiber protein 3, SHIPPO1, Sperm tail protein SHIPPO 1, TISP50
Observed MW:	30 kDa
Uniprot ID :	Q96PU9

Application


Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:500-1:1000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab05973(ODF3 Antibody) at dilution of 1:50



rat testis tissue were subjected to SDS PAGE followed by western blot with FNab05973(ODF3 Antibody) at dilution of 1:600