

anti- OSBPL2 antibody

Product Information

Catalog No.:	FNab06019
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Binds phospholipids; exhibits strong binding to phosphatidic acid and weak binding to phosphatidylinositol 3-phosphate(PubMed:11279184). Binds 25-hydroxycholesterol(PubMed:17428193).

Immunogen information

Immunogen:	oxysterol binding protein-like 2
Synonyms:	KIAA0772, ORP2
Observed MW:	51, 56 kDa
UniprotID :	Q9H1P3

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human ovary tumor using FNab06019(OSBP2 antibody) at dilution of 1:100



A549 cells were subjected to SDS PAGE followed by western blot with FNab06019(OSBP2 Antibody) at dilution of 1:1000