

## anti- OTUD6A antibody

### Product Information

Catalog No.:	FNab06041
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Deubiquitinating enzyme that hydrolyzes 'Lys-27'-, 'Lys-29'-and 'Lys-33'-linked polyubiquitin chains. Also able to hydrolyze 'Lys-11'-linked ubiquitin chains.

### Immunogen information

Immunogen:	OTU domain containing 6A
Synonyms:	DUBA2
Observed MW:	33 kDa
UniprotID :	Q7L8S5

### Application


Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab06041( OTUD6A Antibody) at dilution of 1:1000