

## anti- p120 Catenin antibody

### Product Information

Catalog No.:	FNab06060
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3F8
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	catenin, delta 1
Synonyms:	None
Observed MW:	120 kDa
Uniprot ID :	O60716

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:100-1:400
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human breast ductal carcinoma tissue slide using FNab06060(p120 Catenin antibody) at dilution of 1:200



fetal human brain tissue were subjected to SDS PAGE followed by western blot with FNab06060( p120 Catenin antibody at dilution of 1:1000