

P53 antibody

Product Information


Catalog No.:	FNab06084
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	5D4
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Immunogen:	tumor protein p53
Synonyms:	Antigen NY CO 13, Cellular tumor antigen p53, FLJ92943, LFS1, P53, Phosphoprotein p53, TP53, TRP53, tumor protein p53, Tumor suppressor p53
Observed MW:	53 kDa
Uniprot ID :	P04637

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000-1:10000; IHC: 1:500-1:2000
Image:	



Immunohistochemistry of paraffin-embedded human ovary cancer using FNab06084(P53 antibody) at dilution of 1:500



A431 cells were subjected to SDS PAGE followed by western blot with FNab06084(P53 antibody) at dilution of 1:3000