

anti- PARP15 antibody

Product Information

Catalog No.:	FNab06155
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transcriptional repressor. Has ADP-ribosyltransferase activity.

Immunogen information

Immunogen:	poly(ADP-ribose) polymerase family, member 15
Synonyms:	PARP15, FLJ40196, PARP 15, BAL3, FLJ40597
Observed MW:	50kd
Uniprot ID :	Q460N3

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:1000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



A431 cells were subjected to SDS PAGE followed by western blot with FNab06155(PARP15 antibody) at dilution of 1:500