

## anti- PAX8 antibody

### Product Information

Catalog No.:	FNab06176
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	3F10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Transcription factor for the thyroid-specific expression of the genes exclusively expressed in the thyroid cell type, maintaining the functional differentiation of such cells.

### Immunogen information

Immunogen:	paired box 8
Synonyms:	paired box 8, Paired box protein Pax 8, PAX8
Observed MW:	48 kDa
UniprotID :	Q06710

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:1000-1:4000; IHC: 1:500-1:1000
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab06176(PAX8 Antibody) at dilution of 1:800



SKOV-3 cells were subjected to SDS PAGE followed by western blot with FNab06176(PAX8 Antibody) at dilution of 1:2000