

anti- PDZD3 antibody

Product Information

Catalog No.:	FNab06297
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Acts as a regulatory protein that associates with GUCY2C and negatively modulates its heat-stable enterotoxin-mediated activation. Stimulates SLC9A3 activity in the presence of elevated calcium ions.

Immunogen information

Immunogen:	PDZ domain containing 3
Synonyms:	IKEPP, NHERF4, PDZK2
Observed MW:	54-60KD
UniprotID :	Q86UT5

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Rat kidney tissue were subjected to SDS PAGE followed by western blot with FNab06296(PDZD3 antibody) at dilution of 1:500