

anti- PGLYRP2 antibody

Product Information

Catalog No.:	FNab06356
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May play a scavenger role by digesting biologically active peptidoglycan(PGN) into biologically inactive fragments. Has no direct bacteriolytic activity.

Immunogen information

Immunogen:	peptidoglycan recognition protein 2
Synonyms:	HMFT0141, PGLYRP2, PGLYRPL, PGRP L, PGRPL, tagL, tagL alpha, tagl beta, TAGL like
Observed MW:	68 kDa
Uniprot ID :	Q96PD5

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab06356(PGLYRP2 Antibody) at dilution of 1:50



human plasma tissue were subjected to SDS PAGE followed by western blot with FNab06356(PGLYRP2 Antibody) at dilution of 1:1000