

## anti- PHF5A antibody

### Product Information

Catalog No.:	FNab06387
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Acts as a transcriptional regulator by binding to the GJA1/Cx43 promoter and enhancing its up-regulation by ESR1/ER-alpha. Also involved in pre-mRNA splicing.

### Immunogen information

Immunogen:	PHD finger protein 5A
Synonyms:	None
Observed MW:	12-14 kDa
UniprotID :	Q7RTV0

### Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-PHF5A (IP:FNab06387, 3ug;  
Detection:FNab06387 1:500) with HeLa cells  
lysate 4000ug.



HeLa cells were subjected to SDS PAGE followed  
by western blot with FNab06387(PHF5A  
antibody) at dilution of 1:500