

## anti- PHOX2B antibody

#### **Product Information**

Catalog No.:	FNab06409	
Size:	100µg	
Form:	liquid	
Purification:	Protein A+G purification	
Purity:	$\geq$ 95% as determined by SDS-PAGE	
Host:	Mouse	
Clonality:	monoclonal	
Clone ID:	2B0	
IsoType:	IgG1	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

#### Background

Involved in the development of several major noradrenergic neuron populations, including the locus coeruleus. Transcription factor which could determine a neurotransmitter phenotype in vertebrates. Enhances second-messenger-mediated activation of the dopamine beta-hydrolase and c-fos promoters, and of several enhancers including cAMP-response element and serum-response element.

#### **Immunogen information**

Immunogen:	paired-like homeobox 2b
Synonyms:	NBPhox, Neuroblastoma Phox, paired like homeobox 2b, PHOX2B, PHOX2B homeodomain protein, PMX2B
Observed MW:	35 kDa
UniprotID :	Q99453

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1



# Recommended dilution: WB: 1:1000-1:5000; IHC: 1:20-1:200 Image:





Immunohistochemistry of paraffin-embedded neuroblastoma tissue slide using FNab06409(PHOX2B Antibody) at dilution of 1:50

Neuro-2a cells were subjected to SDS PAGE followed by western blot with FNab06409(PHOX2B Antibody) at dilution of 1:2500

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>