

# anti- PKA C-alpha speicifc antibody

#### **Product Information**

Catalog No.: FNab06478

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

PRKACA, also named as PKACA and PKA C-alpha, belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family and cAMP subfamily. It phosphorylates a large number of substrates in the cytoplasm and the nucleus. This antibody is specific to PRKACA.

## Immunogen information

Immunogen: protein kinase, cAMP-dependent, catalytic, alpha

Synonyms: PKA C alpha, PKACA, PRKACA

Observed MW: 41 kDa Uniprot ID: P17612

## **Application**

Reactivity: Human, Mouse

Tested Application: ELISA, WB, IHC, IP

Recommended dilution: WB: 1:200-1:2000; IP: 1:200-1:2000; IHC: 1:20-1:200

Image:


#### 1

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human testis slide using FNab06478(PRKACA Antibody) at dilution of 1:50



Various lysates were subjected to SDS PAGE followed by western blot with FNab06478(PRKACA antibody) at dilution of 1:500