

anti- PKN0X2 antibody

Product Information

Catalog No.:	FNab06500
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	PBX/knotted 1 homeobox 2
Synonyms:	Homeobox protein PKN0X2, Homeobox protein PREP 2, PBX/knotted 1 homeobox 2, PBX/knotted homeobox 2, PKN0X2, PREP2
Observed MW:	70 kDa,55 kDa
Uniprot ID :	Q96KN3

Application

Reactivity:	Human, Mouse, Rat; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney using FNab06500(PKNOX2 antibody) at dilution of 1:50



RAW264.7 cells were subjected to SDS PAGE followed by western blot with FNab06500(PKNOX2 antibody) at dilution of 1:300