

anti- PLEKHB1 antibody

Product Information

Catalog No.:	FNab06539
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Required for proper localization of retinogeniculate projections but not for eye-specific segregation.

Immunogen information

Immunogen:	pleckstrin homology domain containing, family B(evectins) member 1
Synonyms:	EVT1, KPL1, PHR1, PHRET1
Observed MW:	50 kDa
UniprotID :	Q9UF11

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab06539(PLEKHB1 Antibody) at dilution of 1:300