

anti- POTEA antibody

Product Information

Catalog No.:	FNab06646
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	POTE ankyrin domain family, member A
Synonyms:	A26A1, POTE8
Observed MW:	53 kDa
UniprotID :	Q6S8J7

Application


Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:WB:	1:200-1:1000 ; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab06646(POTEA Antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab06646(POTEA Antibody) at dilution of 1:300