

anti- PPP3R2 antibody

Product Information

Catalog No.:	FNab06732
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Regulatory subunit of calcineurin, a calcium-dependent, calmodulin stimulated protein phosphatase. Confers calcium sensitivity(By similarity).

Immunogen information

Immunogen:	protein phosphatase 3(formerly 2B), regulatory subunit B, beta isoform
Synonyms:	CBLP, PPP3RL
Observed MW:	20 kDa
UniprotID :	Q96LZ3

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab06732(PPP3R2 Antibody) at dilution of 1:200



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab06732(PPP3R2 antibody) at dilution of 1:500