

anti- Protein C inhibitor antibody

Product Information

Catalog No.:	FNab06807
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2H5
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	serpin peptidase inhibitor, clade A(alpha-1 antiproteinase, antitrypsin), member 5
Synonyms:	PAI3, PCI, PLANH3, PROCI, SERPINA5
Observed MW:	46 kDa
Uniprot ID :	P05154

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



human testis tissue were subjected to SDS PAGE followed by western blot with FNab06807(Protein C inhibitor antibody) at dilution of 1:1000