

anti- PRPF19 antibody

Product Information

Catalog No.:	FNab06813
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

PSO4 is the human homolog of yeast Pso4, a gene essential for cell survival and DNA repair (Beck et al., 2008 [PubMed 18263876]).

Immunogen information

Immunogen:	PRP19/PSO4 pre-mRNA processing factor 19 homolog (S. cerevisiae)
Synonyms:	NMP200 ,
Observed MW:	60 kDa
UniprotID :	Q9UMS4

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:100 - 1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded rat brain using FNab06813(PRPF19 Antibody) at dilution of 1:100



SW620 cells were subjected to SDS PAGE followed by western blot with FNab06813(PRPF19 antibody) at dilution of 1:1000