

anti- PSAT1 antibody

Product Information

Catalog No.:	FNab06846
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Catalyzes the reversible conversion of 3-phosphohydroxypyruvate to phosphoserine and of 3-hydroxy-2-oxo-4-phosphonooxybutanoate to phosphohydroxythreonine.

Immunogen information

Immunogen:	phosphoserine aminotransferase 1
Synonyms:	EPIP, Phosphoserine aminotransferase, PSA, PSAT, PSAT1
Observed MW:	35 kDa,40kD
UniprotID :	Q9Y617

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded mouse lung using FNab06846(PSAT1 antibody) at dilution of 1:100



Hela cells were subjected to SDS PAGE followed by western blot with FNab06846(PSAT1 antibody) at dilution of 1:1000