

anti- PSD95-Specific,DLG4 antibody

Product Information

Catalog No.:	FNab06850
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	discs, large homolog 4(Drosophila)
Synonyms:	Disks large homolog 4, DLG4, PSD 95, PSD95, SAP 90, SAP90, Synapse associated protein 90
Observed MW:	90 kDa
UniprotID :	P78352

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:500-1:3000; IP: 1:200-1:1000; IHC: 1:10-1:100
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded rat brain slide using FNab06850(PSD95 Antibody) at dilution of 1:25



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06850(PSD95 Antibody) at dilution of 1:1000