

anti- PTH-Specific antibody

Product Information

Catalog No.:	FNab06923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

PTH belongs to the parathyroid hormone family. PTH elevates calcium level by dissolving the salts in bone and preventing their renal excretion. Defects in PTH are a cause of familial isolated hypoparathyroidism(FIH).

Immunogen information

Immunogen:	parathyroid hormone
Synonyms:	PTH, Parathyrin, Parathormone, PTH1, parathyroid hormone
Observed MW:	18 kDa
Uniprot ID :	P01270

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000-1:4000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Recombinant protein were subjected to SDS PAGE followed by western blot with FNab06923(PTH-Specific antibody) at dilution of 1:2000