

anti- PTPN18 antibody

Product Information

Catalog No.:	FNab06934
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Background

Differentially dephosphorylate autophosphorylated tyrosine kinases which are known to be overexpressed in tumor tissues.

Immunogen information

Immunogen:	protein tyrosine phosphatase, non-receptor type 18(brain-derived)
Synonyms:	BDP1, Brain derived phosphatase, PTP HSCF, PTPN18
Observed MW:	50 kDa
Uniprot ID :	Q99952

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human brain using FNab06934(PTPN18 antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab06934(PTPN18 antibody) at dilution of 1:500