

anti- QRSL1 antibody

Product Information

Catalog No.:	FNab06985
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Allows the formation of correctly charged Gln-tRNA(Gln) through the transamidation of misacylatedGlu-tRNA(Gln) in the mitochondria. The reaction takes place in the presence of glutamine and ATP through an activated gamma-phospho-Glu-tRNA(Gln).

Immunogen information

Immunogen:	glutaminyl-tRNA synthase(glutamine-hydrolyzing)-like 1
Synonyms:	None
Observed MW:	57 kDa
UniprotID :	Q9H0R6

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



human placenta tissue were subjected to SDS
PAGE followed by western blot with
FNab06985(QRSL1 antibody) at dilution of 1:300