

anti- RAB11-Specific antibody

Product Information

Catalog No.: FNab06994

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: RAB11A, member RAS oncogene family

Synonyms: None
Observed MW: 24 kDa
UniprotID: P62491

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IF, IP

Recommended dilution: WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





IP Result of anti-RAB11A-Specific (IP:FNab06994, 3ug; Detection:FNab06994 1:700) with SH-SY5Y cells lysate 3000ug.



MCF-7 cells were subjected to SDS PAGE followed by western blot with 20229-1-AP (RAB11A-Specific antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com