

anti- RAB20 antibody

Product Information

Catalog No.:	FNab07003
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays a role in apical endocytosis/recycling. Plays a role in the maturation and acidification of phagosomes that engulf pathogens, such as *S.aureus* and *M.tuberculosis*. Plays a role in the fusion of phagosomes with lysosomes.

Immunogen information

Immunogen:	RAB20, member RAS oncogene family
Synonyms:	None
Observed MW:	26-29 kDa
UniprotID :	Q9NX57

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-RAB20 (IP:FNab07003, 4ug;
Detection:FNab07003 1:1000) with SW 1990 cells
lysate 2400ug.



SW 1990 cells were subjected to SDS PAGE
followed by western blot with FNab07003(RAB20
antibody) at dilution of 1:1000