

anti- RAB3IP antibody

Product Information

Catalog No.:	FNab07033
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

RAB3IP (RAB3A interacting protein) is a protein-coding gene. GO annotations related to this gene include Rabguanyl-nucleotide exchange factor activity. An important paralog of this gene is RAB3IL1.

Immunogen information

Immunogen:	RAB3A interacting protein (rabin3)
Synonyms:	FLJ14660, FLJ22548, Rab 3A interacting protein, Rab3A interacting protein, RAB3IP, Rabin 3, RABIN3, Rabin8, SSX2 interacting protein
Observed MW:	53 kDa
UniprotID :	Q96QF0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human stomach using FNab07033(RAB3IP antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07033(RAB3IP antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com