

anti- RAB5B antibody

Product Information

Catalog No.:	FNab07040	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)	

Background

Protein transport. Probably involved in vesicular traffic(By similarity).

Immunogen information

Immunogen:	RAB5B, member RAS oncogene family
Synonyms:	RAB5B, Ras related protein Rab 5B
Observed MW:	27 kDa
Uniprot ID :	P61020


Application

Reactivity:	Human, Mouse		
Tested Application:	ELISA, WB		
Recommended dilution: WB: 1:200-1:2000			
Image:			

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab07040(RAB5B antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com