

anti- RAI14 antibody

Product Information

Catalog No.:	FNab07091
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	retinoic acid induced 14
Synonyms:	KIAA1334, NORPEG
Observed MW:	110 kDa
UniprotID :	Q9P0K7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

IP Result of anti-RAI14 (IP:FNab07091, 4ug;
Detection:FNab07091 1:1000) with mouse testis
tissue lysate 4000ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab07091(RAI14 antibody)
at dilution of 1:400

