

anti- RBAP48 antibody

Product Information

Catalog No.:	FNab07149
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	retinoblastoma binding protein 4
Synonyms:	Chromatin assembly factor 1 subunit C,CAF-I 48 kDa subunit,Nucleosome-remodeling factor subunit RBAP48,RBBP-4,Retinoblastoma-binding protein p48
Observed MW:	48-50 kDa
Uniprot ID :	Q09028

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IP, IF
Recommended dilution:	WB: 1:500-1:5000; IP:1:500-1:5000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889


www.fn-test.com



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab07149(RBBP4 antibody) at dilution of 1:100



IP Result of anti-RBAP48 (IP:FNab07149, 4ug; Detection:FNab07149 1:1000) with HeLa cells lysate 2000ug.



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07149(RBBP4 antibody) at dilution of 1:500