

anti- RBAP48 antibody

Product Information

Catalog No.:	FNab07150
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4G10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	retinoblastoma binding protein 4
Synonyms:	Chromatin assembly factor 1 subunit C,CAF-I 48 kDa subunit,Nucleosome-remodeling factor subunit RBAP48,RBBP-4,Retinoblastoma-binding protein p48
Observed MW:	40 kDa, 48-53 kDa
Uniprot ID :	Q09028

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab07150(RBBP4 antibody) at dilution of 1:200



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07150(RBBP4 antibody) at dilution of 1:1000