

anti- RBM19 antibody

Product Information

Catalog No.:	FNab07168
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays a role in embryo pre-implantation development.

Immunogen information

Immunogen:	RNA binding motif protein 19
Synonyms:	DKFZp586F1023, KIAA0682, RBM19, RNA binding motif protein 19
Observed MW:	120kd
Uniprot ID :	Q9Y4C8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human colon using FNab07168(RBM19 antibody) at dilution of 1:50



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab07168(RBM19 antibody) at dilution of 1:300