

anti- RBP4 antibody

Product Information

Catalog No.:	FNab07192
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1C4
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Delivers retinol from the liver stores to the peripheral tissues. In plasma, the RBP-retinol complex interacts with transthyretin, this prevents its loss by filtration through the kidney glomeruli.

Immunogen information

Immunogen:	retinol binding protein 4, plasma
Synonyms:	Plasma retinol binding protein, PRBP, RBP, RBP4, Retinol binding protein 4
Observed MW:	23 kDa
UniprotID :	P02753

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:1000-1:8000; IHC: 1:20-1:200


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human liver cancer slide using FNab07192(RBP4 Antibody) at dilution of 1:50



human blood tissue were subjected to SDS PAGE followed by western blot with FNab07192(RBP4 antibody) at dilution of 1:4000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com