

anti- RCAS1 antibody

Product Information

Catalog No.:	FNab07204
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9E2
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

Background

None

Immunogen information

Immunogen:	estrogen receptor binding site associated, antigen, 9	
Synonyms:	EB9, EBAG9, PDAF, RCAS1	
Observed MW:	28 kDa	
Uniprot ID :	O00559	

Application

Reactivity:	Human		
Tested Application:	ELISA, IHC, WB		
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200			
Image:			


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889


www.fn-test.com





Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab07204(EBAG9 Antibody) at dilution of 1:50 Heat mediated antigen retrieved with Citric acid buffer, pH6.0.

HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07204(EBAG9 Antibody) at dilution of 1:1000



2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com